

Product datasheet for **TS407900P5**

ZNF96 (ZSCAN12) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZSCAN12 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207900
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001039643</u> , <u>NP_001034732</u>
Synonyms:	dj29K1.2; KIAA0426; ZFP96; ZFP96, ZNF305, ZNF29K1, ZSCAN12, KIAA0426, dj29K1.2; zinc finger and SCAN domain containing 12; zinc finger protein 96; zinc finger protein 305; ZNF29K1; ZNF96; ZNF305
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001034732</u>
Locus ID:	9753
Cytogenetics:	6p22.1
Protein Families:	Transcription Factors



[View online »](#)