

Product datasheet for **TS407901**

PFD6 (PFDN6) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PFDN6, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207901
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PFD6 (PFDN6) Mouse Monoclonal Antibody [Clone ID: OTI4E10] (TA504110)
ACCN:	<u>NM_014260</u> , <u>NP_055075</u>
Synonyms:	H2-KE2; HKE2; KE-2; PFD6
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_055075</u>
Locus ID:	10471
Cytogenetics:	6p21.32
Protein Families:	Stem cell - Pluripotency



[View online »](#)