

Product datasheet for **TS408043P5**

PSCD4 (CYTH4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYTH4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208043
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PSCD4 (CYTH4) Mouse Monoclonal Antibody [Clone ID: OT14A5] (TA811273)
ACCN:	<u>NM_013385</u> , <u>NP_037517</u>
Synonyms:	CYT4; cytohesin-4; DJ63G5.1; PSCD4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_037517</u>
Locus ID:	27128
Cytogenetics:	22q13.1



[View online »](#)