

## Product datasheet for **TS408103P5**

### SLAP2 (SLA2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLA2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208103
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SLAP2 (SLA2) Mouse Monoclonal Antibody [Clone ID: OTI1F7] (TA505141)
ACCN:	<u><a href="#">NM_032214</a></u> , <u><a href="#">NP_115590</a></u>
Synonyms:	C20orf156; MARS; SLAP-2; SLAP2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_115590</a></u>
Locus ID:	84174
Cytogenetics:	20q11.23
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)