

## Product datasheet for **TS408406**

### DVL1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DVL1, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208406
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_182779</a></u> , <u><a href="#">NP_877580</a></u>
Synonyms:	dishevelled, dsh homolog 1 (Drosophila); dishevelled 1; dishevelled 1 (homologous to Drosophila dsh); DVL; DVL, MGC54245; MGC54245; OTTHUMP00000003104
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_877580</a></u>
Locus ID:	1855
Cytogenetics:	1p36.33
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS
Protein Pathways:	Basal cell carcinoma, Colorectal cancer, Melanogenesis, Notch signaling pathway, Pathways in cancer, Wnt signaling pathway



[View online »](#)