

## Product datasheet for **TS408422P5**

### LRRC28 CytoSection

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | CytoSections   |
| Description:                          | Transient overexpression of LRRC28 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack  |
| Species:                              | Human  |
| Expression Host:                      | HEK293T  |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC208422   |
| Tag:                                  | C-MYC/DDK  |
| Detection Antibodies:                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)  |
| ACCN:                                 | <u><a href="#">NM_144598</a></u> , <u><a href="#">NP_653199</a></u>  |
| Storage:                              | Room Temperature   |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.  |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.  |
| RefSeq:                               | <u><a href="#">NP_653199</a></u>   |
| Locus ID:                             | 123355   |
| Cytogenetics:                         | 15q26.3  |



[View online »](#)