

Product datasheet for **TS408434P5**

ACBD3 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACBD3 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208434
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ACBD3 Mouse Monoclonal Antibody [Clone ID: OTI5F9] (TA504863)
ACCN:	NM_022735 , NP_073572
Synonyms:	GCP60; GOCAP1; GOLPH1; PAP7
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_073572
Locus ID:	64746
Cytogenetics:	1q42.12
Protein Families:	Druggable Genome



[View online »](#)