

Product datasheet for **TS408469P5**

Upstream Binding Protein 1 (UBP1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UBP1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208469
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Upstream Binding Protein 1 (UBP1) Mouse Monoclonal Antibody [Clone ID: OTI9G7] (TA808851)
ACCN:	NM_014517 , NP_055332
Synonyms:	LBP-1a; LBP-1B; LBP1A; LBP1B
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_055332
Locus ID:	7342
Cytogenetics:	3p22.3
Protein Families:	Transcription Factors



[View online »](#)