

Product datasheet for **TS408586P5**

LAT2 (SLC7A8) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLC7A8, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208586
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	LAT2 (SLC7A8) Mouse Monoclonal Antibody [Clone ID: OTI5A9] (TA500503)
ACCN:	NM_012244 , NP_036376
Synonyms:	LAT2; LPI-PC1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_036376
Locus ID:	23428
Cytogenetics:	14q11.2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)