

Product datasheet for TS408622P5

SRC CytoSection

Product data:

Product Type: CytoSections **Description:** Transient overexpression of SRC, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack Species: Human **Expression Host:** HEK293T Expression cDNA Clone TrueORF Clone RC208622 or AA Sequence: C-MYC/DDK Tag: **Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) ACCN: NM 005417, NP 005408 Synonyms: ASV; c-SRC; p60-Src; SRC1; THC6 Storage: Room Temperature Slides are guaranteed for a year from the date of receipt if proper storage instructions were Stability: followed. **Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. Note: This product is for research use only and is not approved for use in humans or in clinical diagnosis. NP 005408 RefSeq: Locus ID: 6714 Cytogenetics: 20q11.23 **Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase, Stem cell relevant signaling -JAK/STAT signaling pathway **Protein Pathways:** Adherens junction, Endocytosis, Epithelial cell signaling in Helicobacter pylori infection, ErbB signaling pathway, Focal adhesion, Gap junction, GnRH signaling pathway, Tight junction, VEGF signaling pathway



View online 🤉

1/1

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn