

## Product datasheet for **TS408655**

### H3FD (HIST1H3E) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of H3C6 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208655
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_003532</u> , <u>NP_003523</u>
Synonyms:	H3.1; H3/d; H3C1; H3C2; H3C3; H3C4; H3C7; H3C8; H3C10; H3C11; H3C12; H3FD; HIST1H3E
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003523</u>
Locus ID:	8353
Cytogenetics:	6p22.2
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)