

Product datasheet for **TS408666P5**

ROR beta (RORB) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RORB in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208666
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ROR beta (RORB) Mouse Monoclonal Antibody [Clone ID: OT11G1] (TA806996)
ACCN:	NM_006914 , NP_008845
Synonyms:	bA133M9.1; EIG15; NR1F2; ROR-BETA; RZR-BETA; RZRB
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_008845
Locus ID:	6096
Cytogenetics:	9q21.13
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



[View online »](#)