

Product datasheet for **TS408667P5**

UNG CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UNG, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208667
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UNG Mouse Monoclonal Antibody [Clone ID: OTI2C12] (TA503563)
ACCN:	NM_003362 , NP_003353
Synonyms:	DGU; HIGM4; HIGM5; UDG; UNG1; UNG2; UNG15
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003353
Locus ID:	7374
Cytogenetics:	12q24.11
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Base excision repair, Primary immunodeficiency



[View online »](#)