

Product datasheet for **TS408671P5**

ANGPTL1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANGPTL1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208671
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	ANGPTL1 Mouse Monoclonal Antibody [Clone ID: OTI1B2] (TA805608)
ACCN:	NM_004673 , NP_004664
Synonyms:	ANG3; ANGPT3; AngY; ARP1; dj595C2.2; UNQ162
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004664
Locus ID:	9068
Cytogenetics:	1q25.2
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)