

Product datasheet for **TS408752P5**

BCL10 CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of BCL10 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC208752 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) |
| Target Detection Antibodies: | BCL10 Mouse Monoclonal Antibody [Clone ID: OTI4A8] (TA501237) |
| ACCN: | NM_003921 , NP_003912 |
| Synonyms: | c-E10; CARMEN; CIPER; CLAP; IMD37; mE10 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_003912 |
| Locus ID: | 8915 |
| Cytogenetics: | 1p22.3 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | B cell receptor signaling pathway, T cell receptor signaling pathway |



[View online »](#)