

## Product datasheet for TS408838P5

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **ZAP70 CytoSection**

**Product data:** 

**Product Type:** CytoSections

**Description:** Transient overexpression of ZAP70, transcript variant 1, in HEK293T cells, FFPE control for

IHC, ICC and ISH staining, 25 slides per pack

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC208838

Tag: C-MYC/DDK

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

ACCN: <u>NM 001079</u>, <u>NP 001070</u>

Synonyms: ADMIO2; IMD48; SRK; STCD; STD; TZK; ZAP-70

Storage: Room Temperature

**Stability:** Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

**Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

**Note:** This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

RefSeq: NP 001070

 Locus ID:
 7535

 Cytogenetics:
 2q11.2

**Protein Families:** Druggable Genome, Protein Kinase

Protein Pathways: Natural killer cell mediated cytotoxicity, Primary immunodeficiency, T cell receptor signaling

pathway

