

Product datasheet for **TS408860**

PPA2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPA2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC208860
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PPA2 Mouse Monoclonal Antibody [Clone ID: OTI7C1] (TA808324)
ACCN:	NM_006903 , NP_008834
Synonyms:	HSPC124; SCFAI; SCFI; SID6-306
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_008834
Locus ID:	27068
Cytogenetics:	4q24
Protein Pathways:	Oxidative phosphorylation



[View online »](#)