

## Product datasheet for **TS409080P5**

### CRMP2 (DPYSL2) Cytosection

#### Product data:

Product Type:	Cytosections
Description:	Transient overexpression of DPYSL2 (NM_001386) in HEK293T cells paraffin embedded 4 um sections, controls for ICC/IHC staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209080
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	CRMP2 (DPYSL2) Mouse Monoclonal Antibody [Clone ID: OTI5D7] (TA811601)
ACCN:	<a href="#">NM_001386</a> , <a href="#">NP_001377</a>
Synonyms:	CRMP-2; CRMP2; DHPRP2; DRP-2; DRP2; N2A3; ULIP-2; ULIP2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 4 um sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001377</a>
Locus ID:	1808
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome
Protein Pathways:	Axon guidance



[View online »](#)