

## Product datasheet for **TS409163P5**

### AGO1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AGO1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209163
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	AGO1 Mouse Monoclonal Antibody [Clone ID: OTI4G1] (TA811173)
ACCN:	<a href="#">NM_012199</a> , <a href="#">NP_036331</a>
Synonyms:	EIF2C; EIF2C1; GERP95; hAgo1; Q99
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_036331</a>
Locus ID:	26523
Cytogenetics:	1p34.3
Protein Families:	Druggable Genome



[View online »](#)