

Product datasheet for **TS409167P5**

CNKSR3 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CNKSR3 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209167
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CNKSR3 Mouse Monoclonal Antibody [Clone ID: OTI1D7] (TA504737)
ACCN:	NM_173515 , NP_775786
Synonyms:	CNK3; CNK3/IPCEF1; MAGI1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_775786
Locus ID:	154043
Cytogenetics:	6q25.2
Protein Families:	Druggable Genome



[View online »](#)