

Product datasheet for **TS409219P5**

TRIB1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TRIB1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209219
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	TRIB1 Mouse Monoclonal Antibody [Clone ID: OTI8C8] (TA805274)
ACCN:	NM_025195 , NP_079471
Synonyms:	C8FW; GIG-2; GIG2; SKIP1; TRB-1; TRB1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_079471
Locus ID:	10221
Cytogenetics:	8q24.13
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)