

Product datasheet for **TS409226P5**

ADAMTS4 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ADAMTS4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209226
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	ADAMTS4 Mouse Monoclonal Antibody [Clone ID: OTI6F8] (TA805208)
ACCN:	<u>NM_005099</u> , <u>NP_005090</u>
Synonyms:	ADAMTS-2; ADAMTS-4; ADMP-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_005090</u>
Locus ID:	9507
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane



[View online »](#)