

## Product datasheet for **TS409446**

### LYPLA3 (PLA2G15) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PLA2G15 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209446
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_012320</a> , <a href="#">NP_036452</a>
Synonyms:	ACS; GXVPLA2; LLPL; LPLA2; LYPLA3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_036452</a>
Locus ID:	23659
Cytogenetics:	16q22.1
Protein Families:	Secreted Protein
Protein Pathways:	Glycerophospholipid metabolism, Lysosome



[View online »](#)