

## Product datasheet for **TS409568P5**

### CD73 (NT5E) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NT5E, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209568
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD73 (NT5E) Mouse Monoclonal Antibody [Clone ID: OT11G2] (TA809085)
ACCN:	<u><a href="#">NM_002526</a></u> , <u><a href="#">NP_002517</a></u>
Synonyms:	CALJA; CD73; E5NT; eN; eNT; NT; NT5; NTE
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_002517</a></u>
Locus ID:	4907
Cytogenetics:	6q14.3
Protein Families:	ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism, Purine metabolism, Pyrimidine metabolism



[View online »](#)