

## Product datasheet for **TS409579**

### TTDA (GTF2H5) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GTF2H5 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209579
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_207118</u> , <u>NP_997001</u>
Synonyms:	bA120J8.2; C6orf175; TFB5; TFIIH; TGF2H5; TTD; TTD-A; TTD3; TTDA
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_997001</u>
Locus ID:	404672
Cytogenetics:	6q25.3
Protein Families:	Transcription Factors
Protein Pathways:	Nucleotide excision repair



[View online »](#)