

## Product datasheet for **TS409619**

### ACTH (POMC) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POMC, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209619
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ACTH (POMC) Mouse Monoclonal Antibody [Clone ID: OT11A8] (TA506608)
ACCN:	<a href="#">NM_000939</a> , <a href="#">NP_000930</a>
Synonyms:	ACTH; CLIP; LPH; MSH; NPP; OBAIRH; POC
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_000930</a>
Locus ID:	5443
Cytogenetics:	2p23.3
Protein Families:	Druggable Genome
Protein Pathways:	Adipocytokine signaling pathway, Melanogenesis



[View online »](#)