

Product datasheet for **TS409683P5**

Grainyhead like protein 1 homolog (GRHL1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GRHL1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209683
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_014552</u> , <u>NP_055367</u>
Synonyms:	grainyhead-like 1 (Drosophila); LBP-32; LBP32; LBP protein 32; leader-binding protein 32; mammalian grainyhead; MGR; TFCP2L2; transcription factor CP2-like 2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_055367</u>
Locus ID:	29841
Cytogenetics:	2p25.1
Protein Families:	Transcription Factors, Transmembrane



[View online »](#)