

Product datasheet for **TS409701P5**

Influenza Virus NS1A Binding Protein (IVNS1ABP) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IVNS1ABP in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209701
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_006469 , NP_006460
Synonyms:	ARA3; FLARA3; HSPC068; IMD70; KLHL39; ND1; NS-1; NS1-BP; NS1BP
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006460
Locus ID:	10625
Cytogenetics:	1q25.3
Protein Families:	Transcription Factors



[View online »](#)