

## Product datasheet for **TS409846**

### FOXO3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FOXO3, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209846
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FOXO3 Mouse Monoclonal Antibody [Clone ID: OTI5F11] (TA809449)
ACCN:	<a href="#">NM_001455</a> , <a href="#">NP_001446</a>
Synonyms:	AF6q21; FKHL1; FKHL1P2; FOXO2; FOXO3A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001446</a>
Locus ID:	2309
Cytogenetics:	6q21
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Chemokine signaling pathway, Endometrial cancer, Neurotrophin signaling pathway, Non-small cell lung cancer



[View online »](#)