

## Product datasheet for **TS409922P5**

### Ribosomal protein L26 (RPL26) CytoSection

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | CytoSections   |
| Description:                          | Transient overexpression of RPL26 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack   |
| Species:                              | Human  |
| Expression Host:                      | HEK293T  |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC209922   |
| Tag:                                  | C-MYC/DDK  |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)   |
| ACCN:                                 | <u><a href="#">NM_000987</a></u> , <u><a href="#">NP_000978</a></u>  |
| Synonyms:                             | DBA11; L26   |
| Storage:                              | Room Temperature   |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.  |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.  |
| RefSeq:                               | <u><a href="#">NP_000978</a></u>   |
| Locus ID:                             | 6154   |
| Cytogenetics:                         | 17p13.1  |
| Protein Pathways:                     | Ribosome   |



[View online »](#)