

## **Product datasheet for TS410142P5**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Serotonin N acetyltransferase (AANAT) CytoSection

**Product data:** 

**Product Type:** CytoSections

**Description:** Transient overexpression of AANAT, transcript variant 2, in HEK293T cells, FFPE control for

IHC, ICC and ISH staining, 25 slides per pack

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC210142

Tag: C-MYC/DDK

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

Target Detection

Serotonin N acetyltransferase (AANAT) Mouse Monoclonal Antibody [Clone ID: OTI6E1]

Antibodies: (TA810161)

ACCN: <u>NM 001088</u>, <u>NP 001079</u>

Synonyms: DSPS; SNAT

Storage: Room Temperature

**Stability:** Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

**Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

**Note:** This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

**RefSeq:** NP 001079

Locus ID: 15

Cytogenetics: 17q25.1

**Protein Pathways:** Metabolic pathways, Tryptophan metabolism

