

## Product datasheet for **TS410157**

### OR1F1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of OR1F1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210157
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_012360</a></u> , <u><a href="#">NP_036492</a></u>
Synonyms:	OLFMF; OR1F4; OR1F5; OR1F6; OR1F7; OR1F8; OR1F9; OR1F10; OR1F13P; OR3-145; OR16-36; OR16-37; OR16-88; OR16-89; OR16-90; ORL1023
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_036492</a></u>
Locus ID:	4992
Cytogenetics:	16p13.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Olfactory transduction



[View online »](#)