

## Product datasheet for **TS410170**

### Cytochrome P450 3A4 (CYP3A4) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYP3A4, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210170
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_017460</a> , <a href="#">NP_059488</a>
Synonyms:	CP33; CP34; CYP3A; CYP3A3; CYP3A3; CYP3A3; CYP3A3; HLP; NF-25; P450C3; P450PCN1; VDDR3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_059488</a>
Locus ID:	1576
Cytogenetics:	7q22.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, P450, Transmembrane
Protein Pathways:	Drug metabolism - cytochrome P450, Drug metabolism - other enzymes, Linoleic acid metabolism, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Retinol metabolism



[View online »](#)