

Product datasheet for **TS410208P5**

HOXD10 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HOXD10 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210208
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OT14C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	HOXD10 Mouse Monoclonal Antibody [Clone ID: OT11D11] (TA800777)
ACCN:	NM_002148 , NP_002139
Synonyms:	Hox-4.4; HOX4; HOX4D; HOX4E
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002139
Locus ID:	3236
Cytogenetics:	2q31.1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)