

## Product datasheet for **TS410312**

### CNDP1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CNDP1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210312
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CNDP1 Mouse Monoclonal Antibody [Clone ID: OT12F8] (TA500448)
ACCN:	<a href="#">NM_032649</a> , <a href="#">NP_116038</a>
Synonyms:	CN1; CPGL2; HsT2308
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_116038</a>
Locus ID:	84735
Cytogenetics:	18q22.3
Protein Families:	Protease, Secreted Protein
Protein Pathways:	beta-Alanine metabolism, Histidine metabolism, Metabolic pathways



[View online »](#)