

## Product datasheet for **TS410367P5**

### 5HT4 Receptor (HTR4) Cytosection

#### Product data:

Product Type:	Cytosections
Description:	Transient overexpression of HTR4 (NM_000870) in HEK293T cells paraffin embedded 4 um sections, controls for ICC/IHC staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210367
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u><a href="#">NM_000870</a></u> , <u><a href="#">NP_000861</a></u>
Synonyms:	5-HT4; 5-HT4R
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 4 um sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_000861</a></u>
Locus ID:	3360
Cytogenetics:	5q32
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction



[View online »](#)