

Product datasheet for **TS410377P5**

CHST6 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CHST6 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210377
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	CHST6 Mouse Monoclonal Antibody [Clone ID: OTI1D1] (TA800544)
ACCN:	NM_021615 , NP_067628
Synonyms:	C-GlcNAc6ST; gIcNAc6ST-5; gn6st-5; GST4-beta; hCGn6ST; MCDC1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_067628
Locus ID:	4166
Cytogenetics:	16q23.1
Protein Pathways:	Keratan sulfate biosynthesis, Metabolic pathways



[View online »](#)