

Product datasheet for **TS410707P5**

Visfatin (NAMPT) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NAMPT in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210707
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_005746 , NP_005737
Synonyms:	1110035O14Rik; PBEF; PBEF1; VF; VISFATIN
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005737
Locus ID:	10135
Cytogenetics:	7q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Nicotinate and nicotinamide metabolism



[View online »](#)