

Product datasheet for TS410792P5

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PKR (EIF2AK2) CytoSection

Product data:

Product Type: CytoSections

Description: Transient overexpression of EIF2AK2, transcript variant 1, in HEK293T cells, FFPE control for

IHC, ICC and ISH staining, 25 slides per pack

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC210792

Tag: C-MYC/DDK

Detection Antibodies: DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

Target Detection

Antibodies:

PKR (EIF2AK2) Mouse Monoclonal Antibody [Clone ID: OTI3F5] (TA813289)

ACCN: NM 002759, NP 002750

Synonyms: EIF2AK1; LEUDEN; PKR; PPP1R83; PRKR

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

Note: This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

RefSeq: NP 002750

 Locus ID:
 5610

 Cytogenetics:
 2p22.2

Protein Families: Druggable Genome, Protein Kinase, Transcription Factors

