

Product datasheet for **TS410829P5**

MSL2L1 (MSL2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MSL2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210829
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MSL2L1 (MSL2) Mouse Monoclonal Antibody [Clone ID: OTI4G7] (TA811181)
ACCN:	NM_018133 , NP_060603
Synonyms:	MSL-2; MSL2L1; RNF184
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060603
Locus ID:	55167
Cytogenetics:	3q22.3
Protein Families:	Druggable Genome



[View online »](#)