

## Product datasheet for **TS410860P5**

### SAP102 (DLG3) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DLG3, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210860
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SAP102 (DLG3) Mouse Monoclonal Antibody [Clone ID: OT13F2] (TA502190)
ACCN:	<a href="#">NM_021120</a> , <a href="#">NP_066943</a>
Synonyms:	MRX; MRX90; NEDLG; PPP1R82; SAP102; XLMR
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_066943</a>
Locus ID:	1741
Cytogenetics:	Xq13.1



[View online »](#)