

Product datasheet for **TS410961P5**

CD10 (MME) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MME, transcript variant 2a, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210961
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD10 (MME) Mouse Monoclonal Antibody [Clone ID: OT13D11] (TA504507)
ACCN:	NM_007288 , NP_009219
Synonyms:	CALLA; CD10; CMT2T; NEP; SCA43; SFE
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_009219
Locus ID:	4311
Cytogenetics:	3q25.2
Protein Families:	Druggable Genome, Protease, Transmembrane
Protein Pathways:	Alzheimer's disease, Hematopoietic cell lineage, Renin-angiotensin system



[View online »](#)