

Product datasheet for **TS411027P5**

TAS2R48 (TAS2R19) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAS2R19 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211027
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_176888</u> , <u>NP_795369</u>
Synonyms:	MSTP058; T2R19; T2R23; T2R48; TAS2R23; TAS2R48
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_795369</u>
Locus ID:	259294
Cytogenetics:	12p13.2
Protein Families:	Transmembrane
Protein Pathways:	Taste transduction



[View online »](#)