

Product datasheet for **TS411074P5**

SCAND3 (ZBED9) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZBED9 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211074
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_052923 , NP_443155
Synonyms:	Buster4; dj1186N24.3; SCAND3; ZFP38-L; ZNF305P2; ZNF452
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_443155
Locus ID:	114821
Cytogenetics:	6p22.1
Protein Families:	Transcription Factors



[View online »](#)