

## Product datasheet for **TS411151P5**

### HYAL4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HYAL4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211151
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u><a href="#">NM_012269</a></u> , <u><a href="#">NP_036401</a></u>
Synonyms:	CSHY; HYAL-4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_036401</a></u>
Locus ID:	23553
Cytogenetics:	7q31.32
Protein Pathways:	Glycosaminoglycan degradation, Metabolic pathways



[View online »](#)