

Product datasheet for **TS411439P5**

SUCLG1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SUCLG1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211439
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_003849</u> , <u>NP_003840</u>
Synonyms:	GALPHA; MTDPS9; SUCLA1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003840</u>
Locus ID:	8802
Cytogenetics:	2p11.2
Protein Pathways:	Citrate cycle (TCA cycle), Metabolic pathways, Propanoate metabolism



[View online »](#)