

Product datasheet for **TS411653P5**

SPANXA1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SPANXA1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211653
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_013453</u> , <u>NP_038481</u>
Synonyms:	CT11.1; CT11.3; NAP-X; SPAN-X; SPAN-Xa; SPAN-Xb; SPANX; SPANX-A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_038481</u>
Locus ID:	30014
Cytogenetics:	Xq27.2



[View online »](#)