

Product datasheet for **TS411751P5**

RAD50 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RAD50 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211751
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_005732 , NP_005723
Synonyms:	hRad50; NBSLD; RAD502
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005723
Locus ID:	10111
Cytogenetics:	5q31.1
Protein Families:	Druggable Genome
Protein Pathways:	Homologous recombination, Non-homologous end-joining



[View online »](#)