

## Product datasheet for **TS412033**

### CHX10 (VSX2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of VSX2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212033
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_182894</u> , <u>NP_878314</u>
Synonyms:	CHX10; HOX10; MCOP2; MCOPCB3; RET1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_878314</u>
Locus ID:	338917
Cytogenetics:	14q24.3
Protein Families:	Transcription Factors



[View online »](#)