

Product datasheet for **TS412191P5**

Annexin A11 (ANXA11) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANXA11, transcript variant a, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212191
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Annexin A11 (ANXA11) Mouse Monoclonal Antibody [Clone ID: OTI1C6] (TA500950)
ACCN:	NM_001157 , NP_001148
Synonyms:	ALS23; ANX11; CAP-50; CAP50
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001148
Locus ID:	311
Cytogenetics:	10q22.3



[View online »](#)